

# AM3 Integrated Table / Tubestand System



## X-ray System Ideal for Limited Spaces

Integrated table/tubestand saves valuable space and does not require wall or ceiling modifications.

### Features

- X ray equipment compatible with Computed Radiography units (CR) and flat panel detectors (Direct Radiography - DR).
- Unique off-table option provides a full range of off-table procedures including standing knee and ankle exams.
- This ruggedly constructed table will consistently support patient weights up to 350 pounds, making it ideal for busy clinics with demanding caseloads. Heavy-duty, all-steel construction and an extra-long pedestal create a durable table that will last for years.
- Patient positioning is quick and easy on the flat top table. Its smooth surface means no grooves or ridges for clothes to snag or grime to collect.



## Specifications

### Select Generator

- Spectra Series: 20kW, 30kW, 40kW, 50kW

### Integrated Table/Tubestand

- 76" x 30" flat, seamless, stationary table top
- 56" long, enclosed, welded pedestal
- Attached floating tubestand with electric locks
- Angulating platform tube mount with 17" arm
- Angulation dial and operator handgrips
- Electric locks
- 17" grid cabinet with 50" travel
- Grid: 8:1 ratio, 103 LPI
- Grid cabinet to tubestand interlock
- Regular-duty cassette tray
- Operates on 24 VDC
- 300 lb. rated patient load



### Wallstand Package

- Counterbalanced column, floor-to-wall mounted
- Grid cabinet, 17" x 17"
- Grid: 8:1 ratio, 103 LPI
- Regular-duty cassette tray
- Mechanical locks

### X-Ray Tube

- – 2.0mm focal spots, 140,000 heat units
- 0.6-1.5 mm focal spots, 200,000 heat units, specify shoot left or right (with 40 kW generators only)
- 0.6 – 1.2 mm focal spots, 300,000 heat units (with 50 kW generator only)

### Cables

- 20 ft. high voltage cables with federal terminations

### Collimator

- Certified manual collimator with light field and swivel mount

