

# AM4 Stationary Table System



## Versatile X-Ray System Due of its Floor-to-Wall or Ceiling Tubestand

This x-ray equipment has a flat, extra-wide, tabletop that allows easy patient positioning. In addition, the wallstand provides effortless positioning of stationary grid cabinet or optional reciprocating bucky.

### Features

- X ray equipment compatible with Computed Radiography units (CR) and flat panel detectors (Direct Radiography - DR)
- This affordable x-ray system features a sturdy table that supports loads up to 300 pounds. An extra-long pedestal base adds stability and full bucky travel.
- With a detached floor-to-wall (or to ceiling) tubestand, this x-ray equipment is capable of a wide range of radiographic procedures including table-top extremities, skull, and abdominal table work.
- When used with the wallstand, the system capabilities include standing limbs and chest exams.



- This versatile x-ray system may include a stationary fine-line grid or a reciprocating bucky. The cabinet accepts a film cassette tray which accommodates up to a 14" x 17" cassette and 17"x17" flat panel detectors in either direction. Specifications

## Specifications

### Select Generator

- Spectra Series: 20kW, 30kW, 40kW, 50kW

### Stationary Table

- Stationary-top radiographic table, 76" x 30"
- Welded base construction
- Enclosed pedestal
- Grid cabinet with 50" travel
- Grid: 8:1 ratio, 103 LPI
- Regular-duty cassette tray
- 300 lb. rated patient load



### Tubestand

- Floor-to-wall or floor-to-ceiling tubestand with 8 ft. rails
- Electric locks
- Angulation dial and operator handgrips
- Angulating platform tube mounts with 17" arm

### Wallstand Package

- Counterbalanced column, floor-to-wall mounted
- Grid cabinet, 17" x 17"
- Grid: 8:1 ratio, 103 LPI
- Regular-duty cassette tray
- Mechanical locks

### X-Ray Tube

- – 2.0mm focal spots, 140,000 heat units
- 0.6-1.5 mm focal spots, 200,000 heat units, specify shoot left or right (with 40 kW generator only)
- 0.6 – 1.2 mm focal spots, 300,000 heat units (with 50 kW generator only)

### Cables

- 20 ft. high voltage cables with federal terminations

### Collimator

- Certified manual collimator with light field and swivel mount

