

# GP3 Elevating 4 Way Float Top Table System



## World-Class System for Demanding X-Ray Room

The x-ray system GP3 was engineered for performance and long-term reliability. It was designed for efficient workflow in busy environments.

### Features

- Affordability, durability and flexibility
- Excellent x-ray system for hospitals, clinics, orthopedic and general practices.
- Table top supports patients up to 750 lbs. (340 kgs)
- DR, CR or Film solutions

### Specifications

#### High Frequency Generator

- SHF Series: Choice of 20, 32, 40 or 50 kW

#### Elevating, Four Way Float Top Table – 750 lb. Patient Load

- 86.5" x 31.5" float table top
- 43.3" longitudinal and 9.5" lateral table top travel
- Vertical Travel 19.69" x 35.43"
- Collision protection
- Super Speed Bucky



- Heavy duty manual cassette tray
- 103 line, 10:1 ratio, 34-44" focal range grid

**Non-Rotate Tube Stand**

- Floor to wall mounted
- Fits in ceiling heights of 91" to 102"
- Vertical travel of 9.5" to 78"
- Platform tube mount, angulation dial and operator handles
- Magnetic locks control longitudinal, vertical and tube rotation
- 64" Longitudinal Travel
- Rotation option:  $\pm 5$ " Transverse travel, 180° column rotation, vertical travel 14.5" to 76", 80" Longitudinal Travel



**Wall Stand with Grid Cabinet**

- Floor to wall mount with magnetic locks, 84" column
- 17" x 17" grid cabinet with front panel
- Heavy duty manual cassette tray
- 103 line, 10:1 ratio, 40-72" focal range grid

**X-Ray Tube**

- 140,000 heat unit, 1.0-2.0mm focal spots, 90° horn angle
- Other Tubes Optional

**Collimator**

- Certified manual collimator with laser light
- Swivel Mount

**High Voltage Cables**

- 20' (609.6 cm) long

